

Preg Check



Managing the reproductive cycle on your farm

From timely breeding to pregnancy confirmation to successful calving – a critical component on today's dairy farms.

The milk pregnancy test is based on Elisa technology and detects absence of Pregnancy Associated Glycoproteins (PAGs) in the milk sample. It has a 98% level of accuracy.



PregCheck Service includes:-

- Most effective in regular milk recording herds with accurate breeding and service data
- Least invasive pregnancy testing
- Sampling program can be during routine milk recording
- Up to three tests may be defined per service – post service, mid term and pre dry
- First test must be more than 35 days after service and 60 days after calving
- Results may be positive, open and recheck
- Rechecks free of charge
- Access results through Your Herd (www.thecis.co.uk)

Contact us for quotation

01923 695319